



「*Briview*」 Disposable Electronic Flexible Choledochoscope – for ERCP use

Shanghai SeeGen Photoelectric Technology Co., Ltd.



Shanghai SeeGen Photoelectric
Technology Co., Ltd.
Add.: 3 Floor, Building No.1 4299 Jindu
Road, Minhang Dist., Shanghai China
Email: endoexp@163.com

SEEGEN
vision borderless




0123 EN ISO13485:2016

- High resolution image.
- 2 light guide lens make more bright view.



- 1.2mm working channel plus
- 0.6mm water jet channel.



- Easy to insert through duodenoscope working channel.
- Angle-lock provides stable detection procedure.
- Rotate function allows all direction viewing.

Model	CHV-US100J-U
Field of view	110°
Depth of view	5-50mm
Deflection up/down	120°/120°
Insertion dia.	3.0mm
Distal end dia.	3.3mm
Working length	2120mm
Working channel	1.2mm/0.6mm



Optional fixing system to choledochoscope

ROI: Unit 31 Southern Cross Business Park, Boghall Road, Bray, Co Wicklow, A98 PX00
 Tel: +353 1 201 4620
 email: sales@fleetwoodhealthcare.ie

NI: Innis House, Holywood, Belfast, Co. Down, BT18 9HF Tel: +44 28 9099 5160

FLEETWOOD
 HEALTHCARE EXCELLENCE THROUGH INNOVATION

SeeCholERCPO1/20