

BRIVIEW

SINGLE USE CHOLEDOCHOSCOPE

- ▶ Tapered-Shape distal end aids insertion
- ▶ Larger working channel available for added flexibility
- ▶ HD image quality on par with conventional scopes

Universal Processor - compatible with range of disposable scopes



Available in 3 sizes for Surgical & Laparoscopic Applications

Model	CHV 100-E	CHV-100T-E	CHV-US100-H
Field of View	110°	110°	110°
Depth of View	5-50mm	5-50mm	5-50mm
Deflection Up/Down	160°/130°	160°/130°	120°/120°
Insertion Diameter	4.7mm	5.3mm	3.0mm
Distal End Diameter	4.7mm	5.3mm	3.0mm
Working Length	38.5cm	38.5cm	67cm
Working Channel	2.2mm	2.6mm	1.2mm

Brought to you by:

CE 0123 EN ISO13485:2016

SEEGEN 视介光电
vision borderless

Benefits of Briview

Single-Use Choledochoscope

1

Safety - A new sterile scope for every patient, reducing the risk of infection

2

Consistency - Predicable clinical performance everytime - no issue with wear and tear

3

24 x 7 Availability - No need to wait for reprocessing of reusable scopes

4

Cost Transparency - No Service or Leasing Agreements

5

Maintenance Free - **No** degradation, **No** reprocessing and **No** repairs

6

Compatibility - works with most commonly used diagnostic and therapeutic tools

Distributed in Ireland by:

FLEETWOOD
HEALTHCARE
EXCELLENCE THROUGH
INNOVATION

BRVIEWCHOL2020REV1

Fleetwood Healthcare

ROI: Unit 31 Southern Cross Business Park, Boghall Road, Bray, Co. Wicklow, A98 PX00 **Tel:** +353 1 201 4620

NI: Innis House, Holywood, Co. Down, BT18 9HF
Tel: +44 28 9099 5160

e-mail: sales@fleetwoodhealthcare.ie

Website: www.fleetwoodhealthcare.ie