



「*Briview*」 Disposable Electronic Flexible Choledochoscope

Shanghai SeeGen Photoelectric Technology Co., Ltd.


Shanghai SeeGen Photoelectric
Technology Co.,Ltd.
Add.: 3 Floor, Building No.1 4299 Jindu
Road, Minhang Dist., Shanghai China
Email:endoexp@163.com

SEEGEN
vision borderless



0123 EN ISO13485:2016

- Single use avoids cross infection and saves sterilization costs.
- Bigger working channel allows bigger instrument use.
- Taper-shape distal tip makes insertion easier.
- High resolution image makes procedure safer and easier.



Model	CHV-100-E	CHV-100T-E	CHV-US100-H
Field of view	110°	110°	110°
Depth of view	5-50mm	5-50mm	5-50mm
Deflection up/down	210°/130°	210°/130°	120°/120°
Insertion dia.	4.7mm	5.3mm	3.0mm
Distal end dia.	4.7mm	5.3mm	3.0mm
Working length	385mm	385mm	670mm
Working channel	2.2mm	2.6mm	1.2mm



ROI: Unit 31 Southern Cross Business Park, Boghall Road, Bray, Co Wicklow, A98 PX00
 Tel: +353 1 201 4620
 email: sales@fleetwoodhealthcare.ie

NI: Innis House, Holywood, Belfast, Co. Down, BT18 9HF Tel: +44 28 9099 5160

FLEETWOOD
 HEALTHCARE EXCELLENCE THROUGH INNOVATION