


# BRIVIEW

## SINGLE USE CHOLEDOCHOSCOPE FOR ERCP

- HD image quality on par with reusable scopes
- Dual-light guide enhances visualisation
- 1.2mm working channel with 0.6mm water jet channel
- Rotate function & angle lock allows for precise positioning
- Easy to insert through the duodenoscope working channel

Universal Processor - compatible with range of disposable scopes



Model	CHV-US100J-U
Field of View	110°
Depth of View	5-50mm
Deflection Up/Down	120°/120°
Insertion Diameter	3.0mm
Distal End Diameter	3.0mm
Working Length	212cm
Working Channel	1.2mm/0.6mm

Brought to you by:

# Benefits of Briview

Single-Use Choledochoscope for ERCP

1

**Safety** - A new sterile scope for every patient, reducing the risk of infection

2

**Consistency** - Predicable clinical performance everytime - no issue with wear and tear

3

**24 x 7 Availability** - No need to wait for reprocessing of reusable scopes

4

**Cost Transparency** - No Service or Leasing Agreements

5

**Maintenance Free** - **No** degradation, **No** reprocessing and **No** repairs

6

**Compatibility** - works with most commonly used diagnostic and therapeutic tools

Distributed in Ireland by:

**FLEETWOOD**  
HEALTHCARE  
EXCELLENCE THROUGH  
INNOVATION

BRVIEWERCP2020REV1

**Fleetwood Healthcare**

**ROI:** Unit 31 Southern Cross Business Park, Boghall Road, Bray, Co. Wicklow, A98 PX00 **Tel:** +353 1 201 4620

**NI:** Innis House, Holywood, Co. Down, BT18 9HF  
**Tel:** +44 28 9099 5160

**e-mail:** [sales@fleetwoodhealthcare.ie](mailto:sales@fleetwoodhealthcare.ie)

**Website:** [www.fleetwoodhealthcare.ie](http://www.fleetwoodhealthcare.ie)