

Endoscopic
Set

Made in
Germany

New
Device

HaemoCerTM Endoscopic SET

Haemostatic powder for
endoscopic application



Tube length 2400 mm



Outer tube diameter of 2.5 mm



HEPA filter -
High Efficiency Particulate Air




Usable instrument channel:
2.8 mm and above



BioCer

Endoscopic SET

Haemostatic Powder + Endoscopic Applicator



HaemoCer™ Haemostatic Powder

- Absorbable Polysaccharide Haemostat
- Containing no animal nor human components
- Plant-based polymer crosslinking that creates ultra-hydrophilic, biocompatible particles
- Rapidly resorbed by amylase within 48h without leaving any residues behind
- Rapid off-the-shelf deployment, no special storage, defrosting or premixing



HaemoCer™ Endoscopic Applicator

- For upper and lower gastrointestinal bleedings
- Even effective in difficult to access regions
- No blocking in case of contact with moisture
- Tube length 2400 mm
- Outer tube diameter of 2.5 mm
- HEPA filter - High Efficiency Particulate Air
- Usable instrument channel: 2.8 mm and above
- No high pressure gas source needed
- No compressor system needed
- Reduced dust formation
- Precise dosing of the applied powder

Product Overview HaemoCer™ Endoscopic Set

HFS323 HaemoCer™ PLUS 3 g + Endoscopic Applicator 1 pcs



BioCer

Fleetwood Healthcare

Unit 31, Southern Cross Business Park, Boghall Road, Bray,, Co. Wicklow

Tel: (353) 1 2014620 | Website: www.fleetwoodhealthcare.ie | Email: sales@fleetwoodhealthcare.ie



Management
System
EN ISO
13485:2016

www.tuv.com
ID 0000059855



HDD144.02/2020-03