


# UROFINO

## SINGLE USE CYSTOSCOPE

- ▶ Tapered-Shape distal end aids insertion
- ▶ Dual light guide enhances visualisation
- ▶ Larger working channel available for added flexibility
- ▶ HD image quality on par with conventional scopes



Universal Processor - compatible with range of disposable scopes



**Available in 2 sizes with the option of larger working Channel**

Model	UV-100-E	UV-100-T-E
Field of View	110°	110°
Depth of View	5-50mm	5-50mm
Deflection Up/Down	210°/130°	210°/130°
Insertion Diameter	4.7mm	5.3mm
Distal End Diameter	4.7mm	5.3mm
Working Length	38.5cm	38.5cm
Working Channel	2.2mm	2.6mm

Brought to you by:

CE 0123 EN ISO13485:2016

**SEEGEN** 视介光电  
vision borderless

# Benefits of Urofino

Single-Use Cystoscope

1

**Safety** - A new sterile scope for every patient, reducing the risk of infection

2

**Consistency** - Predicable clinical performance everytime - no issue with wear and tear

3

**24 x 7 Availability** - No need to wait for reprocessing of reusable scopes

4

**Cost Transparency** - No Service or Leasing Agreements

5

**Maintenance Free** - **No** degradation, **No** reprocessing and **No** repairs

6

**Compatibility** - works with most commonly used diagnostic and therapeutic tools

Distributed in Ireland by:

**FLEETWOOD**  
HEALTHCARE  
EXCELLENCE THROUGH  
INNOVATION

UROFINOCYSTO2020REV1

**Fleetwood Healthcare**

**ROI:** Unit 31 Southern Cross Business Park, Boghall Road, Bray, Co. Wicklow, A98 PX00 **Tel:** +353 1 201 4620

**NI:** Innis House, Holywood, Co. Down, BT18 9HF  
**Tel:** +44 28 9099 5160

**e-mail:** [sales@fleetwoodhealthcare.ie](mailto:sales@fleetwoodhealthcare.ie)

**Website:** [www.fleetwoodhealthcare.ie](http://www.fleetwoodhealthcare.ie)