


UROFINO

SINGLE USE URETEROSCOPE

- ▶ Tapered-Shape distal end aids insertion
- ▶ Dual light guide enhances visualisation
- ▶ HD image quality on par with conventional scopes
- ▶ Lightweight and easy to control

Universal Processor - compatible with range of disposable scopes



275° bending angle
for improved vision

Model	UV-US100-H
Field of View	110°
Depth of View	5-50mm
Deflection Up/Down	275°/275°
Insertion Diameter	3.0mm
Distal End Diameter	3.0mm
Working Length	75cm
Working Channel	1.2mm

Brought to you by:

CE 0123 EN ISO13485:2016

SEEGEN 视介光电
vision borderless

Benefits of Urofino

Single-Use Ureteroscope

1

Safety - A new sterile scope for every patient, reducing the risk of infection

2

Consistency - Predicatable clinical performance everytime - no issue with wear and tear

3

24 x 7 Availability - No need to wait for reprocessing of reusable scopes

4

Cost Transparency - No Service or Leasing Agreements

5

Maintenance Free - **No** degradation, **No** reprocessing and **No** repairs

6

Compatibility - works with most commonly used diagnostic and therapeutic tools

Distributed in Ireland by:

FLEETWOOD
HEALTHCARE
EXCELLENCE THROUGH
INNOVATION

UROFINO2020 REV1

Fleetwood Healthcare

ROI: Unit 31 Southern Cross Business Park, Boghall Road, Bray,
Co. Wicklow, A98 PX00 **Tel:** +353 1 201 4620

NI: Innis House, Holywood, Co. Down, BT18 9HF
Tel: +44 28 9099 5160

e-mail: sales@fleetwoodhealthcare.ie

Website: www.fleetwoodhealthcare.ie